



Multi-function body pressure dispersion prevents bedsore

The mattress uses "Especially high-density soft polyurethane foam". Its dual structure with moderate resiliency gently supports the body to prevent the bedsores. It connects to the air bathtub with an inward drainage system.

Liberty Himawari (automatic excretion treatment apparatus) is also available with this device.

Air bathtub, easy to set up and put away

The air bathtub can be set up in the bed by just shifting the user's posture without moving them out of bed. Avoiding the accidents, like falling off the bathtub. Its compact storage allows a single caregiver to carry, set up and put away in a moment.

Even a 180cm-tall person can stretch out.

Cover Sheet for water splash prevention

Set up frame poles as a bed guard, with a cover sheet preventing water splash over the walls or floor. The user can bathe at ease within a room.







How to set up the air bathtub

Prepare the air bathtub on the mattress. Care giver can set it up while shifting the user's body, just as changing the bedsheet.







[Specifications of Ajisai's Main Unit]

- Dimensions of the main device: Length: 1,910mm (internal dimension 1,860mm)
 / Width: 800mm (internal dimension 600mm)
 / Height: 300mm
- Weight of the main unit: Approximately 3kg
 Materials: Polyurethane and vinyl chloride